


# Visualize:Vascular Specifications



Salient Imaging provides electronic transfer of the images and the reports via its cloud. Cloud access needs to be installed on one of the PCs in your facility to do this. An authorized Salient Imaging representative will perform the installation for you.

## PC Requirements:

**Windows® 7 Professional or XP professional 32 bit operating system**  
**2GB RAM**  
**10 GB free space on hard drive**  
**Your PC must be networked with Images Source and Destination Addresses**  
**Port 22 needs to be opened in your firewall to send information out**

## Image Requirements:

Images transported to PC via USB (DICOM Part 10) or DICOM network

## Other Information:

**Visualize:Vascular™ accepts this image format:**  
**DICOM multiframe**

**Visualize:Vascular™ generates reports in these formats:**  
**PDF**  
**DICOM PDF**  
**DICOM SC**

**Agency: FDA 510(k)**

## Compatible with the following Ultrasound Equipment:

Ultrasound machine that produces DICOM formatted images.  
Ultrasound machine that can capture at least a 7 sec cine clip, 10 sec cine clip is preferred.  
Ultrasound machine with a pixel size at 80 micron or smaller, note: this excludes the HP5500 and other machines from this era.

|            |  |
|------------|--|
| Esaote     | MyLab® 70, MyLab® 50, MyLab® 40                              |
| GE         | Vivid®Q, Vivid® 7, Logiq® 7, Logiq® 9, Logiqbook®, Logiq® E9 |
| Medison    |  |
| Mindray    |  |
| Philips    | iE 33, HD®11, iU®22, HDI® 5000, HDI® 3500                    |
| Siemens    | S2000®, Sequoia®, Cypress                                    |
| Sonosite   |  |
| Terason    |  |
| Toshiba    | Aplio®, Aplio® XG, Aplio® XU                                 |
| Ultrasonix | SonixTouch®, SonixTablet®                                    |
| Zonare     |  |

© 2013, Salient Imaging, Inc., *Visualize:Vascular* is a trademark of Salient Imaging.. Acrobat is a registered trademark of Adobe. Vivid, Logiq are registered trademarks of GE. Windows and SQL Server are registered trademarks of Microsoft. Sequoia, S2000, Cypress are registered trademarks of Siemens. iU, HD, HDI, and iE are registered trademarks of Phillips. MyLab is a registered trademark of Esaote. Aplio is a registered trademark of Toshiba. SonixTouch and SonixTablet are registered trademarks of Ultrasonix. All rights reserved. Specifications subject to change without notice.

**Salient Imaging**

One Elm Street  
Milford, NH 03055, USA  
603-272-6066  
www.salientimaging.com